



Product Name: DERMOTIC EAR & SKIN SUSPENSION
APVMA Approval No: 53950/112414

Label Name:	DERMOTIC EAR & SKIN SUSPENSION
Signal Headings:	PRESCRIPTION ANIMAL REMEDY KEEP OUT OF REACH OF CHILDREN FOR ANIMAL TREATMENT ONLY
Constituent Statements:	23 mg/ml MICONAZOLE NITRATE 5 mg/ml PREDNISOLONE ACETATE 0.696 mg/ml POLYMYXIN B SULFATE
Claims:	For the treatment of otitis externa and skin infections caused by fungi, yeasts, Gram-negative, and Gram-positive bacteria in horses, dogs and cats.
Net Contents:	20mL 30mL 40mL
Directions for Use:	Shake well before use. Dermotic Ear & Skin Suspension is not active against ear mites.
Restrains:	Not applicable
Contraindications:	Not applicable
Precautions:	Not applicable
Side Effects:	Not applicable
Dosage and Administration:	Otitis externa: Instil 3 - 5 drops into the ear twice daily. Gently massage the ear to ensure even distribution of the product. Skin infections: Apply to the infected area and rub in. After clinical signs have completely disappeared, continue treatment without interruption for several days. In certain obstinate

cases, it may be necessary to continue treatment for 2 - 3 weeks. At the beginning of treatment, cut the hair surrounding or covering lesions.

General Directions:	Not applicable
---------------------	----------------

Withholding Periods:	MEAT WITHHOLDING PERIOD (HORSES): Nil
----------------------	---------------------------------------

Trade Advice:	Not applicable
---------------	----------------

Safety Directions:	Not applicable
--------------------	----------------

First Aid Instructions:	If poisoning occurs, contact a doctor or Poisons Information Centre. Phone Australia 131126.
-------------------------	--

First Aid Warnings:	Not applicable
---------------------	----------------

Additional User Safety:	Not applicable
-------------------------	----------------

Environmental Statements:	Not applicable
---------------------------	----------------

Disposal:	Dispose of empty container by wrapping with paper and putting in the garbage.
-----------	---

Storage:	Store below 25 °C (air conditioning).
----------	---------------------------------------